

permobil

THE POWER OF MOBILITY™



DYNAMIC BY DESIGN

Featuring a rear-wheel drive power base and two different speed packages of up to 6.5 MPH, the C350 Corpus 3G can easily keep pace with your busy lifestyle. Plus, the Corpus 3G seating system is designed for comfort beyond comparison, something you need in today's fast-paced world. The seating system follows the contours of the body, offers infinite adjustments for a custom fit and reclines 150°. The armrests feature built-in channels for cable routing, which greatly reduces the need for messy wire ties; while the legrests are curved for better support. Blending form and function, the C350 Corpus 3G is the ideal blueprint for mobility.

C350
CORPUS 3G





POWER TILT

50° of tilt with two range selections for valuable repositioning.

RIDE COMFORT

4-wheel adjustable independent suspension for the ultimate ride quality.

POWER RECLINE

Up to 150° of recline to find "just the right level" of comfort.

INDEPENDENT ARMRESTS

Adjustable armrest height, angle and pitch as well as a convenient flip-up feature.

AVAILABLE SHROUD COLORS



ONYX BLACK



CRYSTAL BLUE



PLATINUM SILVER



EXTREME YELLOW



VINTAGE WINE

SPECIFICATIONS

Maximum User Weight	300 lbs	136 kg	Seat to Floor Height (Elevator)	17.5"	445 mm
Maximum Speed Options	5.5 MPH, 6.5 MPH	8 km/h, 10.5 km/h	(Seat Tube)	17.5", 18.25", 19"	445, 465, 485 mm
Driving Range ¹	22 Miles	35 km	Tilt Options (Power)	0° – 50°	
Base Width	24.5"	620 mm	(Power)	-5° – 45°	
Base Length	36"	915 mm	Recline Options (Power)	85° – 175°	
Minimum Turning Radius	26"	660 mm	(Manual)	85° – 120°	
Ground Clearance	3"	80 mm	Legrest Elevation Options (Power)	85° – 170°	
Optional Seat Elevator	8"	200 mm	(Manual)	85° – 170°	
Battery Type	Group 34		Armrest Pad Lengths	10", 13", 16"	260, 335, 410 mm
Drive Electronics	90A / 120A R-net		Seat Widths	17", 19", 21", 23"	420, 470, 520, 570 mm
Maximum Slope ²	10° Up / 10° Down		Seat Depths	14" – 22"	370 – 570 mm
			Backrest Heights	20", 23" – 28"	500, 570 – 710 mm

¹ ACTUAL DRIVING RANGE WILL VARY BASED ON DRIVING AND BATTERY CONDITIONS.
² REFER TO OWNER MANUAL FOR SAFE DRIVING PRACTICES